

# Protocol for Talos L120C imaging

## Location:

DE-782 (Talos room)

## General rules:

- Do not use the Talos unless you are accompanied by an EM core staff member or have been cleared as an independent user
- Use tools only for their intended purposes
- Remove grids from desk promptly
- Respect your reservation time
- Report any problems to Caleigh ASAP

## Sign-up protocol:

- Use iLab (base rate, \$60/hr) to reserve for the entire time you will be using the Talos
- Book supported use unless you are **entirely** comfortable with the setup for screening/collection
- Staff services can be booked under “Services” through iLab
- Door access will be granted after training

## Training plan:

- Please contact staff at [emsr@fredhutch.org](mailto:emsr@fredhutch.org) to schedule training
  - Two 1-hour Talos training sessions
  - If returning >30 days after last Talos use, please contact staff and booked “Supported Use” time

## Training objectives:

- Successfully insert the room temperature sample holder without disrupting the column vacuum
- Find eucentric height, focus, and set appropriate imaging conditions
- Convert files to tifs and transfer to lab storage

## User provided materials:

- Sample grids
- Storage space for images

## Shared resources tools list: if anything is missing/out of the ordinary, please contact Caleigh ASAP

- Sharp tweezers
- Room temperature sample holder and pin tool
- Data stored for 30-days after collection

### Startup:

- Turn on the [Filament](#) – allow approx 13min to saturate (status: beam emission stable)
- Sign in at Kiosk and click Start to begin recording filament time
- Insert specimen rod (**\*EM staff will do this unless user is trained\***)
- Set apertures
  - for biological samples [C2 Aperture](#) – 100, [Objective aperture](#) – 70, [Spot Size](#) 6
- Open [Column Valves \(Setup tab\)](#)
- Insert FLU screen ([Screen Up/Down](#))
- Find beam and move to area away from ROI to reduce beam damage (screen on left monitor)
- Center beam in LM and M/SA.  
Under [Direct Alignments](#) tab, select [Beam Shift](#), center beam with [MF X / Y](#), press [Done](#) to save for session  
Quick alignment can be done with [Roller Ball](#)
- Determine Eucentric (Z) Height for each sample using a recognizable feature
  - Mag @ 13,500x, [Wobbler](#) on (flapout on [Search](#) tab), set recognizable feature to swing less than 500nm using [Z +/-](#) button. Easiest to do on FLU screen (Mesh grids close to zero, slot grids +80-85)

### Imaging:

- Joystick moves sample stage, +/- settings to increase or decrease stage movement speed
- Focus image
  - True focus can be determined using the [FFT](#) function (either on the FLU screen or with camera). Click [Reset Defocus](#) to save Eucentric.
  - Focus by setting the [focus step](#) (bottom collar on focus knob) to 4 or below. The higher the focus step number, the coarser the focus step. Focus the image by turning the [focus knob](#) (top) until FFT image is diffuse (no Thon rings). Most biologicals read best a little under focused. Underfocus by turning [focus](#) counter clockwise until image appears sharp but not grainy (**can you recommend a value here?**).
- Camera
  - Open [Camera](#) tab, [Insert](#) camera, retract FLU screen (if not already retracted), select [Search](#). Once ROI is determined and in focus ([FFT](#)), select [Preview](#) for acquisition resolution image (can fine tune focus here), select [Acquire](#) to save image (see Saving Images set up below).

### Saving Images to a Folder:

- Create a new folder on the D drive, under TIA data: (e.g. yourHutchNetID\_today's date)

- To setup Autosave **Insert** the Camera, select **Preview** and then deselect **Preview** to activate **Autosave**.
- Go to the **Autosave** icon on TIA software screen (or under **File** select **autosave**) and click until you get the dialog box, browse to find your folder in D:TIA Data  
To save numbered images, make sure **Save Sequentially** is checked (bottom left)
- Open folder (can drag box to move to left monitor).
- Select the file name box (.emi ), enter the title for image series (e.g. SM\_5432\_2).
- **Acquire** image and click **Save** in the dialog box. Check to see that it goes into selected folder and/or subfolder and is titled correctly. Continue to acquire images until you are done with that sample.
- Repeat for each new sample.

### Batch Converting Images:

- Blank the beam (Camera tab, **Blank**) while you are doing the batch convert
- Select **Folder Export** from left sidebar on TIA screen, open **Settings**
- To save converted images to the same folder, the “Source” and “Target” line names must be the same:
- In TIA Data select “yourFolder” (or subfolder), activate and copy the path line from the top left. Paste into both Source and Target
- Select Image Format as “PC tiff with scale marker (full Res)” to include scale bar
- Say “OK” (data box disappears unless there is an error with the path).
- Click “Export” to start conversion.
- Once everything is converted to tiffs, “Close All (pages icon with x)”.
- You will get a dialog box – Say “No to all”
- Close all documents and discard changes = “Yes”

### How to Separate TIFFs:

- Open Folder and/or subfolder
- View
- Group by
- Type
- Control, Shift – click on first and last image and move TIFFs into a TIFF folder that you have made in the subfolder or into homelinks

### Insert the new sample:

- Click off autosave icon before you focus for another image
- Focus, and **Preview**
- Turn **Preview** off

- Click on [Autosave](#)
- [File](#) – [Autosave](#) (click more than one time) and it will come up with the folder you were in. Make a new subfolder and click on it.
- Go to file name and put in everything before .emi. Ex. SM\_5432\_3
- Don't press "Save" in the dialog box until you have acquired an image.
- Acquire image
- Check to see if the images are in the subfolder.

#### End of Session:

- Close [Column Valves \(Setup tab\)](#)
- Reset [Holder](#) (flapout on [Search](#) tab), zeroing X,Y,Z
- Turn off [AutoSave](#) (!!)
- Retract Camera (deselect [Insert Camera](#))
- Insert FLU screen
- Retract [Obj Ap](#)
- [FLU](#) screen on "Natural"
- Turn off [Filament](#)
- Remove sample – done by EM staff unless trained
- End Kiosk session and Log Out