

Experimental Histopathology

Research AdministrationSeattle, WA ● 501(c)(3) Nonprofit

Fred Hutch's Shared Resources are catalysts for lifesaving discoveries. This uniquely centralized program of 15 specialized core facilities and scientific services drives advances by integrating dedicated experts and cutting-edge technologies across the entire research pipeline, from basic science to clinical trial.

Hardware:

- Leica LMD6
- Laser type: Diode pumped, solid state
- Wavelength: 355nm
 Avg. pulse energy: 70µJ
 Repetition Rate: 80Hz

Camera

- Leica DFC7000 T
- Brightfield and fluorescence (DAPI, FITC, up to Cy5.5)
- More information about Leica DFC7000 T

Camera Sensor

- Sony® ICX674AQG CCD, EXview HAD II technology, quad-tap
- 2.8 megapixels
- Exposure time 4 μsec 200 sec

Scope/Laser

• Detailed technical specifications

Capabilities

- Single cell isolation
- Precise analysis in routine and custom antibody validations
- Isolation of cells and tissues from both FFPE and frozen (OCT) sections
- Isolation of cells and tissues on H&E, chromogenic DAB IHC, and immunofluorescent slides
- Any species



Experimental Histopathology Amanda Koehne, D.V.M., Ph.D., DACVP Director 206.667.6971 akoehne@fredhutch.org



